






## Product data sheet REF 64.54.35

11.12.2025

General data			
article number	64.54.35	diamond burr 95mm 3.5mm	
Unit	piece	Customs code (HS code)	90189084
GMDN	46392	GTIN	04046826109493
Risk category	2a (EU MDD)	UMDNS-Code	15-883
Classification	6 (EU MDD)	Base UDI-DI	
Net weight (kg)	0,004	Operating instructions	
Latex free	yes	Preparation instructions	MPK02
			
Medicon media	450.00.46 (catalogue sinuscopy )		
Complete picture		Detail images	
 			
size	95mm 3.5mm	Coating	Artificially produced diamond coating on working ends of instruments
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part	1.4112 (X90CrMoV18)		