



Product data sheet REF 66.29.37

11.12.2025

General data			
article number	66.29.37	cupped fcps 13cm 90° cvd upw 2mm horiz LL	
Unit	piece	Customs code (HS code)	90189084
GMDN	62467	GTIN	04046826111489
Risk category	1r (EU MDR)	UMDNS-Code	15-800
Classification	6 (EU MDR)	Base UDI-DI	4046826;B;R003;06;ai;5X
Net weight (kg)	0,150	Operating instructions	GATDR00X
Latex free	yes	Preparation instructions	MPK08-1



Medicon media	450.00.37 (catalogue ENT) 450.00.46 (catalogue sinuscopy)
---------------	--

Complete picture

Detail images



size	13cm 90° 2mm Ø 2.7mm	Coating	
------	----------------------	---------	--

Storage conditions

No special requirements

Materials

	Material number (DIN designation)
Main part	1.4021 (X20Cr13)
Screw	1.4305 (X8CrNiS18-9)
Handle	1.4021 (X20Cr13)
Others	1.4301 (X5CrNi18-10)