








Product data sheet REF 66.29.57

11.12.2025

General data			
article number	66.29.57	cupped fcps 13cm 70° cvd upw 3mm vert LL	
Unit	piece	Customs code (HS code)	90189084
GMDN	62467	GTIN	04046826111540
Risk category	1r (EU MDR)	UMDNS-Code	15-800
Classification	6 (EU MDR)	Base UDI-DI	4046826;B;R003;06;ai;5X
Net weight (kg)	0,150	Operating instructions	GATDR00X
Latex free	yes	Preparation instructions	MPK08-1
 			
Medicon media	450.00.37 (catalogue ENT) 450.00.46 (catalogue sinuscopy)		
Complete picture		Detail images	
  			
size	13cm 70° 3mm Ø 2.7mm	Coating	
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part	1.4021 (X20Cr13)		
Screw	1.4305 (X8CrNiS18-9)		
Handle	1.4021 (X20Cr13)		
Others	1.4301 (X5CrNi18-10)		