







## Product data sheet REF 66.33.21

11.12.2025

General data			
article number	66.33.21	punch 16cm Ø 3.5mm 60° ang circular cut	
Unit	piece	Customs code (HS code)	90189084
GMDN	35541	GTIN	04046826112240
Risk category	1r (EU MDR)	UMDNS-Code	13-230
Classification	6 (EU MDR)	Base UDI-DI	4046826;B;R003;06;ai;5X
Net weight (kg)	0,150	Operating instructions	GATDR00X
Latex free	yes	Preparation instructions	MPK08-1
 			
Medicon media	450.00.37 (catalogue ENT ) 450.00.46 (catalogue sinuscopy )		
Complete picture		Detail images	
 			
size	16cm Ø 3.5mm60°	Coating	
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
<b>Main part</b>	1.4021 (X20Cr13)		
<b>Screw</b>	1.4305 (X8CrNiS18-9)		
<b>Handle</b>	1.4021 (X20Cr13)		
<b>Others</b>	1.4301 (X5CrNi18-10)		