








## Product data sheet REF 66.38.21

11.12.2025

| General data  |   |                               |                         |
|---|---|-------------------------------|-------------------------|
| article number  | 66.38.21  | Hajek punch 13cm 90° up 3.5mm |                         |
| Unit  | piece   | Customs code (HS code)        | 90189084                |
| GMDN  | 35541   | GTIN                          | 04046826112325          |
| Risk category   | 1r (EU MDR)   | UMDNS-Code                    | 13-230                  |
| Classification  | 6 (EU MDR)  | Base UDI-DI                   | 4046826;B;R003;06;ai;5X |
| Net weight (kg)   | 0,152   | Operating instructions        | GATDR00X                |
| Latex free  | yes   | Preparation instructions      | MPK23-1                 |
|     |   |                               |                         |
| Medicon media   | 450.00.37 (catalogue ENT )<br>450.00.46 (catalogue sinuscopy )<br>450.00.62 (catalogue Medicon instruments Oral and Maxillofacial Surgery)<br>451.02.51 (ENT basic sets Dr M W Winter Instr sets f ENT surgery) |                               |                         |
| Complete picture  |   | Detail images                 |                         |
|    |   |                               |                         |
| size  | 13cm 90° 3.5mm  | Coating                       |                         |
| Storage conditions  |   |                               |                         |
| No special requirements   |   |                               |                         |
| Materials   |   |                               |                         |
|   | Material number (DIN designation)   |                               |                         |
| Main part   | 1.4021 (X20Cr13)  |                               |                         |
| Screw   | 1.4305 (X8CrNiS18-9)  |                               |                         |
| Handle  | 1.4021 (X20Cr13)  |                               |                         |
| Others  | 1.4310 (X10CrNi18-8)  |                               |                         |