







Product data sheet REF 66.38.22

11.12.2025

General data			
article number	66.38.22	Hajek punch 13cm 90° dw 3.5mm	
Unit	piece	Customs code (HS code)	90189084
GMDN	35541	GTIN	04046826112332
Risk category	1r (EU MDR)	UMDNS-Code	13-230
Classification	6 (EU MDR)	Base UDI-DI	4046826;B;R003;06;ai;5X
Net weight (kg)	0,146	Operating instructions	GATDR00X
Latex free	yes	Preparation instructions	MPK23-1
 			
Medicon media	450.00.37 (catalogue ENT) 450.00.46 (catalogue sinuscopy) 450.00.62 (catalogue Medicon instruments Oral and Maxillofacial Surgery)		
Complete picture		Detail images	
 			
size	13cm 90° 3.5mm	Coating	
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part	1.4021 (X20Cr13)		
Screw	1.4305 (X8CrNiS18-9)		
Handle	1.4021 (X20Cr13)		
Others	1.4310 (X10CrNi18-8)		