








## Product data sheet REF 66.38.41

11.12.2025

General data			
article number	66.38.41	McKenty sphenoid bone punch thrcut 2mm	
Unit	piece	Customs code (HS code)	90189084
GMDN	35541	GTIN	04046826112394
Risk category	1r (EU MDR)	UMDNS-Code	13-230
Classification	6 (EU MDR)	Base UDI-DI	4046826;B;R003;06;ai;5X
Net weight (kg)	0,095	Operating instructions	GATDR00X
Latex free	yes	Preparation instructions	MPK23-1
 			
Medicon media	450.00.37 (catalogue ENT ) 450.00.46 (catalogue sinuscopy )		
Complete picture		Detail images	
  			
size	2mm 17cm	Coating	
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
<b>Main part</b>	1.4021 (X20Cr13)		
<b>Screw</b>	1.4305 (X8CrNiS18-9)		
<b>Handle</b>	1.4021 (X20Cr13)		
<b>Others</b>	1.4310 (X10CrNi18-8)		