








## Product data sheet REF 66.54.10

11.12.2025

General data			
article number	66.54.10	Killian septum elevator d/e 18cm	
Unit	piece	Customs code (HS code)	90189084
GMDN	35047	GTIN	04046826112745
Risk category	1r (EU MDR)	UMDNS-Code	11-506
Classification	6 (EU MDR)	Base UDI-DI	4046826;B;R004;05;aj;6A
Net weight (kg)	0,024	Operating instructions	GATDR00X
Latex free	yes	Preparation instructions	MPK03-1
 			
Medicon media	450.00.37 (catalogue ENT )		
Complete picture		Detail images	
  			
size	18cm	Coating	
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part	1.4021 (X20Cr13)		