








Product data sheet REF 66.54.50

11.12.2025

General data			
article number	66.54.50	Pennington septum elevator d/e 20cm	
Unit	piece	Customs code (HS code)	90189084
GMDN	35047	GTIN	04046826112790
Risk category	1r (EU MDR)	UMDNS-Code	11-506
Classification	6 (EU MDR)	Base UDI-DI	4046826;B;R004;05;aj;6A
Net weight (kg)	0,038	Operating instructions	GATDR00X
Latex free	yes	Preparation instructions	MPK03-1
 			
Medicon media	450.00.37 (catalogue ENT) 450.00.60 (catalogue plastic surgery)		
Complete picture		Detail images	
  			
size	20cm	Coating	
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part	1.4021 (X20Cr13)		