








## Product data sheet REF 66.59.70

11.12.2025

General data			
article number	66.59.70	Sewall raspatory 6mm 20cm	
Unit	piece	Customs code (HS code)	90189084
GMDN	35047	GTIN	04046826113148
Risk category	1r (EU MDR)	UMDNS-Code	11-506
Classification	6 (EU MDR)	Base UDI-DI	4046826;B;R003;03;ai;4W
Net weight (kg)	0,025	Operating instructions	GATDR00X
Latex free	yes	Preparation instructions	MPK20-1
 			
Medicon media	450.00.37 (catalogue ENT ) 450.00.46 (catalogue sinuscopy )		
Complete picture		Detail images	
  			
size	6mm 20cm	Coating	
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part	1.4021 (X20Cr13)		
Others	1.4305 (X8CrNiS18-9)		