



Product data sheet REF 66.69.15

11.12.2025

General data			
article number	66.69.15	flex nose pharyngoscope Ø 3.8mm 300mm 8000 pixles	
Unit	piece	Customs code (HS code)	90189084
GMDN	36709	GTIN	04046826113568
Risk category	1 (EU MDR)	UMDNS-Code	13-011
Classification	5 (EU MDR)	Base UDI-DI	4046826;N004;01;XD
Net weight (kg)	3,300	Operating instructions	GATDN00X
Latex free	yes	Preparation instructions	
 			
Medicon media	450.00.46 (catalogue sinuscopy)		
Complete picture		Detail images	
size	Ø 3.8mm 300mm 9800 Pixel	Coating	
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part	KEINE ANGABEN (Werkstoffangabe nicht vorhanden bzw. Stückliste)		