



Product data sheet REF 66.70.40

11.12.2025

| General data | | | |
|---|--|----------------------------|----------------|
| article number | 66.70.40 | endoscope Ø 4.0mm 175mm 0° | |
| Unit | piece | Customs code (HS code) | 90189020 |
| GMDN | 37180 | GTIN | 04046826113773 |
| Risk category | 2a (EU MDD) | UMDNS-Code | 13-396 |
| Classification | 5, 6 (EU MDD) | Base UDI-DI | |
| Net weight (kg) | 0,160 | Operating instructions | G776000 |
| Latex free | yes | Preparation instructions | |
|  | | | |
| Medicon media | 450.00.37 (catalogue ENT) 450.00.46 (catalogue sinuscopy) | | |
| Complete picture | | Detail images | |
|  | | | |
| size | Ø 4.0mm 175mm 0° | Coating | |
| Storage conditions | | | |
| No special requirements | | | |
| Materials | | | |
| | Material number (DIN designation) | | |
| Main part | 1.4301 (X5CrNi18-10) | | |
| Screw | 1.4305 (X8CrNiS18-9) | | |
| Others | PEEK (Polyetheretherketon) | | |