








Product data sheet REF 67.15.35

11.12.2025

General data			
article number	67.15.35	Yankauer septum elevator 18cm 6mm	
Unit	piece	Customs code (HS code)	90189084
GMDN	35047	GTIN	04046826114398
Risk category	1r (EU MDR)	UMDNS-Code	11-506
Classification	6 (EU MDR)	Base UDI-DI	4046826;B;R004;05;aj;6A
Net weight (kg)	0,033	Operating instructions	GATDR00X
Latex free	yes	Preparation instructions	MPK03-1
 			
Medicon media	450.00.37 (catalogue ENT) 450.00.60 (catalogue plastic surgery)		
Complete picture		Detail images	
  			
size	18cm 6mm	Coating	
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part Handle	1.4021 (X20Cr13) 1.4305 (X8CrNiS18-9)		