








Product data sheet REF 67.16.55

11.12.2025

| General data | | | |
|---|--|---------------------------------------|-------------------------|
| article number | 67.16.55 | Cottle McKenty raspatory d/e 4mm 22cm | |
| Unit | piece | Customs code (HS code) | 90189084 |
| GMDN | 35559 | GTIN | 04046826114602 |
| Risk category | 1r (EU MDR) | UMDNS-Code | 11-506 |
| Classification | 6 (EU MDR) | Base UDI-DI | 4046826;B;R003;03;ai;4W |
| Net weight (kg) | 0,038 | Operating instructions | GATDR00X |
| Latex free | yes | Preparation instructions | MPK03-1 |
|   | | | |
| Medicon media | 450.00.03 (general catalogue the art of surgery) 450.00.37 (catalogue ENT) 450.00.46 (catalogue sinuscopy) 450.00.60 (catalogue plastic surgery) | | |
| Complete picture | | Detail images | |
|    | | | |
| size | 4mm 22cm | Coating | |
| Storage conditions | | | |
| No special requirements | | | |
| Materials | | | |
| | Material number (DIN designation) | | |
| Main part | 1.4021 (X20Cr13) | | |