








## Product data sheet REF 67.20.05

11.12.2025

General data			
article number	67.20.05	Kilner retr sh 10cm	
Unit	piece	Customs code (HS code)	90189084
GMDN	45918	GTIN	04046826114701
Risk category	1r (EU MDR)	UMDNS-Code	12-033
Classification	6 (EU MDR)	Base UDI-DI	4046826;B;R004;03;aj;5L
Net weight (kg)	0,015	Operating instructions	GATDR00X
Latex free	yes	Preparation instructions	MPK01-1
 			
Medicon media	450.00.03 (general catalogue the art of surgery) 450.00.37 (catalogue ENT ) 450.00.60 (catalogue plastic surgery )		
Complete picture		Detail images	
  			
size	10cm	Coating	
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part	1.4021 (X20Cr13)		