








Product data sheet REF 67.20.21

11.12.2025

| General data | | | |
|---|--|---------------------------------|-------------------------|
| article number | 67.20.21 | Cottle retr 2t 14cm ball tip/sh | |
| Unit | piece | Customs code (HS code) | 90189084 |
| GMDN | 45918 | GTIN | 04046826114749 |
| Risk category | 1r (EU MDR) | UMDNS-Code | 12-033 |
| Classification | 6 (EU MDR) | Base UDI-DI | 4046826;B;R004;03;aj;5L |
| Net weight (kg) | 0,020 | Operating instructions | GATDR00X |
| Latex free | yes | Preparation instructions | MPK03-1 |
|   | | | |
| Medicon media | 450.00.03 (general catalogue the art of surgery) 450.00.37 (catalogue ENT) 450.00.60 (catalogue plastic surgery) | | |
| Complete picture | | Detail images | |
|    | | | |
| size | 14cm | Coating | |
| Storage conditions | | | |
| No special requirements | | | |
| Materials | | | |
| | Material number (DIN designation) | | |
| Main part | 1.4021 (X20Cr13) | | |