







## Product data sheet REF 67.20.25

11.12.2025

General data			
article number	67.20.25	Fomon retr 2t ball tips 17cm	
Unit	piece	Customs code (HS code)	90189084
GMDN	45918	GTIN	04046826114770
Risk category	1r (EU MDR)	UMDNS-Code	12-033
Classification	6 (EU MDR)	Base UDI-DI	4046826;B;R004;03;aj;5L
Net weight (kg)	0,021	Operating instructions	GATDR00X
Latex free	yes	Preparation instructions	MPK03-1
 			
Medicon media	450.00.37 (catalogue ENT ) 450.00.60 (catalogue plastic surgery )		
Complete picture		Detail images	
 			
size	17cm	Coating	
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part Handle	1.4021 (X20Cr13) 1.4024 (X15Cr13)		