








Product data sheet REF 67.60.63

11.12.2025

General data			
article number	67.60.63	Fomon HM rasp fig 7/8 20cm	
Unit	piece	Customs code (HS code)	90189084
GMDN	35559	GTIN	04046826116316
Risk category	1r (EU MDR)	UMDNS-Code	13-290
Classification	6 (EU MDR)	Base UDI-DI	4046826;B;R003;03;ai;4W
Net weight (kg)	0,055	Operating instructions	GATDR00X
Latex free	yes	Preparation instructions	MPK12-1
 			
Medicon media	450.00.03 (general catalogue the art of surgery) 450.00.37 (catalogue ENT)		
Complete picture		Detail images	
  			
size	20cm	Coating	Artificially created galvanic gold layer (max 0.04mm) on parts of the surface of instruments
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part	1.4021 (X20Cr13)		
Others	K40 (WcCo)		