



Product data sheet REF 67.75.07

11.12.2025

| General data | | | |
|---|--|---|-------------------------|
| article number | 67.75.07 | Entner ligature carrier f mini lift 7cm | |
| Unit | piece | Customs code (HS code) | 90189084 |
| GMDN | 32864 | GTIN | 04046826116965 |
| Risk category | 1r (EU MDR) | UMDNS-Code | 12-729 |
| Classification | 6 (EU MDR) | Base UDI-DI | 4046826;B;R001;03;ag;3J |
| Net weight (kg) | 0,030 | Operating instructions | GATDR00X |
| Latex free | yes | Preparation instructions | |
|   | | | |
| Medicon media | 450.00.60 (catalogue plastic surgery) | | |
| Complete picture | | Detail images | |
|  | | | |
| size | 7cm | Coating | |
| Storage conditions | | | |
| No special requirements | | | |
| Materials | | | |
| | Material number (DIN designation) | | |
| Main part Handle | 1.4021 (X20Cr13) 1.4301 (X5CrNi18-10) | | |