







Product data sheet REF 68.02.78

11.12.2025

| General data | | | |
|--|--|--|-------------------------|
| article number | 68.02.78 | intra oral retr 28cm w fibre optic illum | |
| Unit | piece | Customs code (HS code) | 90189084 |
| GMDN | 45918 | GTIN | 04046826118082 |
| Risk category | 1r (EU MDR) | UMDNS-Code | 12-028 |
| Classification | 6 (EU MDR) | Base UDI-DI | 4046826;B;R004;03;aj;5L |
| Net weight (kg) | 0,150 | Operating instructions | GATDR00X |
| Latex free | yes | Preparation instructions | |
|   | | | |
| Medicon media | 450.00.03 (general catalogue the art of surgery) 450.00.62 (catalogue Medicon instruments Oral and Maxillofacial Surgery) | | |
| Complete picture | | Detail images | |
|   | | | |
| size | 28cm | Coating | |
| Storage conditions | | | |
| No special requirements | | | |
| Materials | | | |
| | Material number (DIN designation) | | |
| Main part | 1.4021 (X20Cr13) | | |
| Handle | 1.4301 (X5CrNi18-10) | | |
| Others | GLASFASER (Dünne Faser die aus Glas besteht) | | |