







Product data sheet REF 68.03.29

11.12.2025

| General data | | | |
|--|---|---------------------------------------|-------------------------|
| article number | 68.03.29 | Obwegeser maxilla mobilizer 23cm 12mm | |
| Unit | piece | Customs code (HS code) | 90189084 |
| GMDN | 12846 | GTIN | 04046826118211 |
| Risk category | 1r (EU MDR) | UMDNS-Code | 12-846 |
| Classification | 6 (EU MDR) | Base UDI-DI | 4046826;B;R003;01;ai;48 |
| Net weight (kg) | 0,124 | Operating instructions | GATDR00X |
| Latex free | yes | Preparation instructions | MPK20-1 |
|   | | | |
| Medicon media | 450.00.03 (general catalogue the art of surgery) 450.00.60 (catalogue plastic surgery) 450.00.62 (catalogue Medicon instruments Oral and Maxillofacial Surgery) 451.09.02 (brochure bimaxillary surgery) | | |
| Complete picture | | Detail images | |
|   | | | |
| size | 23cm 12mm | Coating | |
| Storage conditions | | | |
| No special requirements | | | |
| Materials | | | |
| | Material number (DIN designation) | | |
| Main part | 1.4021 (X20Cr13) | | |