







## Product data sheet REF 68.04.73

11.12.2025

General data			
article number	68.04.73	Obwegeser HM wire twister 12.5cm	
Unit	piece	Customs code (HS code)	90189084
GMDN	32886	GTIN	04046826118655
Risk category	1r (EU MDR)	UMDNS-Code	15-221
Classification	6 (EU MDR)	Base UDI-DI	4046826;B;R002;02;ah;3V
Net weight (kg)	0,030	Operating instructions	GATDR00X
Latex free	yes	Preparation instructions	
 			
Medicon media	450.00.03 (general catalogue the art of surgery) 450.00.60 (catalogue plastic surgery ) 450.00.62 (catalogue Medicon instruments Oral and Maxillofacial Surgery)		
Complete picture		Detail images	
 			
size	12.5cm	Coating	Artificially created galvanic gold layer (max 0.04mm) on parts of the surface of instruments
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part	1.4021 (X20Cr13)		
Others	K40 (WcCo)		