



Product data sheet REF 68.23.12A

25.08.2026

General data			
article number	68.23.12A	Ti cross slot screw Ø2.4 L=12mm self tapping with locking	
Unit	box/10	Customs code (HS code)	90211090
GMDN	46638	GTIN	04046826537289
Risk category	(EU MDD)	UMDNS-Code	
Classification	NA (EU MDD)	Base UDI-DI	4046826;B001;01;P8
Net weight (kg)	0,000	Operating instructions	G689480
Latex free	yes	Preparation instructions	MPK06
			
Medicon media			
Complete picture		Detail images	
size	2.4 12mm	Coating	
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		