









Product data sheet REF 68.24.40

11.12.2025

General data			
article number	68.24.40	hand piece ang f Kirschner wires	
Unit	piece	Customs code (HS code)	90189084
GMDN	44491	GTIN	04046826119287
Risk category	N.A. (EU MDD)	UMDNS-Code	
Classification	NA (EU MDD)	Base UDI-DI	
Net weight (kg)	0,239	Operating instructions	
Latex free		Preparation instructions	KEIN MP
			
Medicon media			
Complete picture		Detail images	
Storage conditions			
 60 °C	 90 %	 1060 hPa	
	 700 hPa		
Materials			
	Material number (DIN designation)		
Main part	KEINE ANGABEN (Werkstoffangabe nicht vorhanden bzw. Stückliste)		