







Product data sheet REF 68.51.02

11.12.2025

General data			
article number	68.51.02	Tessier maxillary mobilizer rt	
Unit	piece	Customs code (HS code)	90189084
GMDN	33432	GTIN	04046826120016
Risk category	1r (EU MDR)	UMDNS-Code	12-028
Classification	6 (EU MDR)	Base UDI-DI	4046826;B;R004;03;aj;5L
Net weight (kg)	0,095	Operating instructions	GATDR00X
Latex free	yes	Preparation instructions	MPK12-1
 			
Medicon media	450.00.03 (general catalogue the art of surgery) 450.00.62 (catalogue Medicon instruments Oral and Maxillofacial Surgery)		
Complete picture		Detail images	
 			
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part	1.4021 (X20Cr13)		
Screw	1.4305 (X8CrNiS18-9)		