







Product data sheet REF 68.70.40

11.12.2025

General data			
article number	68.70.40	Ti Mini L plate 100° lf p/2	
Unit	box/2	Customs code (HS code)	90211090
GMDN	46642	GTIN	04046826120238
Risk category	2b (EU MDD)	UMDNS-Code	13-050
Classification	8 (EU MDD)	Base UDI-DI	
Net weight (kg)	0,010	Operating instructions	G687004
Latex free	yes	Preparation instructions	MPK06
  			
Medicon media			
Complete picture		Detail images	
			
size	100°	Coating	
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part	3.7035 (ASTM_Grade_2)		