







## Product data sheet REF 68.70.44

11.12.2025

General data			
article number	68.70.44	Ti Mini L plate 100° If int sect 12mm p/2	
Unit	box/2	Customs code (HS code)	90211090
GMDN	46642	GTIN	04046826120276
Risk category	2b (EU MDD)	UMDNS-Code	13-050
Classification	8 (EU MDD)	Base UDI-DI	
Net weight (kg)	0,011	Operating instructions	G687004
Latex free	yes	Preparation instructions	MPK06
  			
Medicon media			
Complete picture		Detail images	
			
size	100° 12mm	Coating	
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part	3.7035 (ASTM_Grade_2)		