







Product data sheet REF 68.73.37A

11.12.2025

General data			
article number	68.73.37A	Ti Mini L plate rt int sect 15mm Low Prof p/2	
Unit	box/2	Customs code (HS code)	90211090
GMDN	46642	GTIN	04046826342128
Risk category	2b (EU MDD)	UMDNS-Code	13-050
Classification	8 (EU MDD)	Base UDI-DI	
Net weight (kg)	0,010	Operating instructions	G687004
Latex free	yes	Preparation instructions	MPK06
  			
Medicon media			
Complete picture		Detail images	
			
size	15mm	Coating	
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part	3.7035 (ASTM_Grade_2)		