







Product data sheet REF 68.74.14

11.12.2025

| General data | | | |
|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-------------------------------------------|----------------------------------------------|----------------|
| article number | 68.74.14 | Ti Mini plate 4 holes short int sect 8mm p/2 | |
| Unit | box/2 | Customs code (HS code) | 90211090 |
| GMDN | 46642 | GTIN | 04046826121204 |
| Risk category | 2b (EU MDD) | UMDNS-Code | 13-050 |
| Classification | 8 (EU MDD) | Base UDI-DI | |
| Net weight (kg) | 0,009 | Operating instructions | G687004 |
| Latex free | yes | Preparation instructions | MPK06 |
|    | | | |
| Medicon media | 451.09.02 (brochure bimaxillary surgery) | | |
| Complete picture | | Detail images | |
|  | | | |
| size | 8mm | Coating | |
| Storage conditions | | | |
| No special requirements | | | |
| Materials | | | |
| | Material number (DIN designation) | | |
| Main part | 3.7035 (ASTM_Grade_2) | | |