







## Product data sheet REF 68.74.70A

11.12.2025

General data			
article number	68.74.70A	Ti Mini L plate 110° rt+lf int sect 10mm Rigid p/2	
Unit	box/2	Customs code (HS code)	90211090
GMDN	46642	GTIN	04046826417147
Risk category	2b (EU MDD)	UMDNS-Code	13-050
Classification	8 (EU MDD)	Base UDI-DI	
Net weight (kg)	0,000	Operating instructions	G687004
Latex free	yes	Preparation instructions	MPK06
  			
Medicon media			
Complete picture		Detail images	
			
size	110° 10mm	Coating	
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part	3.7035 (ASTM_Grade_2)		