







Product data sheet REF 68.75.15A

11.12.2025

General data			
article number	68.75.15A	Ti Mini slot screw 2.0x15mm grey self-tapp p/10	
Unit	box/10	Customs code (HS code)	90211090
GMDN	61837	GTIN	04046826342845
Risk category	2b (EU MDD)	UMDNS-Code	16-101
Classification	8 (EU MDD)	Base UDI-DI	
Net weight (kg)	0,010	Operating instructions	G687004
Latex free	yes	Preparation instructions	MPK06
  			
Medicon media			
Complete picture		Detail images	
			
size	2.0x15mm	Coating	
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part	3.7165 (ASTM_Grade_5_TiAl6V4_(ELI))		