







Product data sheet REF 68.75.42A

11.12.2025

General data			
article number	68.75.42A	Ti Mini Emerg Easy Grip screw 2.3x5.5mm red self-tapp p/10	
Unit	box/10	Customs code (HS code)	90211090
GMDN	61837	GTIN	04046826343125
Risk category	2b (EU MDD)	UMDNS-Code	16-101
Classification	8 (EU MDD)	Base UDI-DI	
Net weight (kg)	0,009	Operating instructions	G687004
Latex free	yes	Preparation instructions	MPK06
  			
Medicon media	451.08.13 (brochure intraoral distraction devices f upper a lower jaws)		
Complete picture		Detail images	
			
size	2.3x5.5mm	Coating	
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part	3.7165 (ASTM_Grade_5_TiAl6V4_(ELI))		