





## Product data sheet REF 68.76.33

11.12.2025

General data			
article number	68.76.33	twist drill Ø 1.2x8mm 50mm Stryker	
Unit	piece	Customs code (HS code)	90184910
GMDN	10520	GTIN	04046826122041
Risk category	2a (EU MDD)	UMDNS-Code	11-341
Classification	6 (EU MDD)	Base UDI-DI	
Net weight (kg)	0,002	Operating instructions	G930100
Latex free	yes	Preparation instructions	MPK02
			
Medicon media	451.08.13 (brochure intraoral distraction devices f upper a lower jaws)		
Complete picture		Detail images	
			
size	Ø 1.2x8mm 50mm	Coating	
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part	1.4112 (X90CrMoV18)		