



Product data sheet REF 68.76.37

25.08.2026

General data			
article number	68.76.37	twist drill Ø 1.2x6.5mm 101.5mm Stryker	
Unit	piece	Customs code (HS code)	90184910
GMDN	32390	GTIN	04046826542184
Risk category	1r (EU MDR)	UMDNS-Code	11-331
Classification	6 (EU MDR)	Base UDI-DI	4046826;B;R003;05;ai;5L
Net weight (kg)	0,007	Operating instructions	GATDR00X
Latex free	yes	Preparation instructions	MPK01-1
			
Medicon media			
Complete picture		Detail images	
size	Ø 1.2x6.5mm 101.5mm	Coating	
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part	1.4112 (X90CrMoV18)		