





Product data sheet REF 68.76.45

11.12.2025

General data			
article number	68.76.45	twist drill Ø 1.45x25mm 115mm Stryker	
Unit	piece	Customs code (HS code)	90184910
GMDN	10520	GTIN	04046826401436
Risk category	2a (EU MDD)	UMDNS-Code	11-341
Classification	6 (EU MDD)	Base UDI-DI	
Net weight (kg)	0,000	Operating instructions	G930100
Latex free	yes	Preparation instructions	MPK02
			
Medicon media			
Complete picture		Detail images	
			
size	Ø 1.45x25mm 115mm	Coating	Color thinner
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part	1.4034 (X46Cr13)		