








Product data sheet REF 68.76.57

11.12.2025

General data			
article number	68.76.57	depth gauge f drill-holes 20cm f syst 2.4	
Unit	piece	Customs code (HS code)	90189084
GMDN	32845	GTIN	04046826122263
Risk category	1r (EU MDR)	UMDNS-Code	15-237
Classification	6 (EU MDR)	Base UDI-DI	4046826;B;R005;02;ak;5X
Net weight (kg)	0,050	Operating instructions	GATDR00X
Latex free	yes	Preparation instructions	MPK12-1
 			
Medicon media			
Complete picture		Detail images	
  			
size	20cm	Coating	
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part	1.4301 (X5CrNi18-10)		
Others	1.4310 (X10CrNi18-8)		