





Product data sheet REF 68.76.73

11.12.2025

General data			
article number	68.76.73	twist drill Ø 1.45x11mm 44mm Stryker	
Unit	piece	Customs code (HS code)	90184910
GMDN	10520	GTIN	04046826336707
Risk category	2a (EU MDD)	UMDNS-Code	11-341
Classification	6 (EU MDD)	Base UDI-DI	
Net weight (kg)	0,000	Operating instructions	G930100
Latex free	yes	Preparation instructions	MPK02
			
Medicon media	451.08.13 (brochure intraoral distraction devices f upper a lower jaws) 451.09.02 (brochure bimaxillary surgery) 453.06.00_DEU (Hybrid Plate System german Brochure) 453.06.00_ENG (Hybrid Plate System - Reconstruction of the mandible with a microvascular fibular bone graft ENG Brochure) 453.06.00_ESP (Hybrid Plate System - Reconstruction of the mandible with a microvascular fibular bone graft ESP Brochure)		
Complete picture		Detail images	
			
size	Ø 1.45x11mm 44mm	Coating	
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part	1.4112 (X90CrMoV18)		