





Product data sheet REF 68.78.01

11.12.2025

| General data | | | |
|--|-----------------------------------|------------------------------------|----------------|
| article number | 68.78.01 | twist drill Ø 1.45x4mm 25mm dental | |
| Unit | piece | Customs code (HS code) | 90184910 |
| GMDN | 10520 | GTIN | 04046826336752 |
| Risk category | 2a (EU MDD) | UMDNS-Code | 11-341 |
| Classification | 6 (EU MDD) | Base UDI-DI | |
| Net weight (kg) | 0,000 | Operating instructions | G930100 |
| Latex free | yes | Preparation instructions | MPK02 |
|  | | | |
| Medicon media | 453.12.00 (Epiplating) | | |
| Complete picture | | Detail images | |
|  | | | |
| size | Ø 1.45x4mm 25mm | Coating | |
| Storage conditions | | | |
| No special requirements | | | |
| Materials | | | |
| | Material number (DIN designation) | | |
| Main part | 1.4112 (X90CrMoV18) | | |