



## Product data sheet REF 68.78.03

11.12.2025

General data			
article number	68.78.03	twist drill Ø 1.45x5mm 34mm dental	
Unit	piece	Customs code (HS code)	90184910
GMDN	10520	GTIN	04046826336769
Risk category	2a (EU MDD)	UMDNS-Code	11-341
Classification	6 (EU MDD)	Base UDI-DI	
Net weight (kg)	0,000	Operating instructions	G930100
Latex free	yes	Preparation instructions	MPK02
			
Medicon media	453.12.00 (Epiplating )		
Complete picture		Detail images	
			
size	Ø 1.45x5mm 34mm	Coating	
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part	1.4112 (X90CrMoV18)		