




Product data sheet REF 68.78.93

11.12.2025

| General data | | | |
|---|-----------------------------------|--------------------------------|----------------|
| article number | 68.78.93 | Twist drill Ø1.45x18mm 27.06mm | |
| Unit | piece | Customs code (HS code) | 90184910 |
| GMDN | 10520 | GTIN | 04046826554538 |
| Risk category | 2a (EU MDD) | UMDNS-Code | 11-341 |
| Classification | 6 (EU MDD) | Base UDI-DI | |
| Net weight (kg) | 0,000 | Operating instructions | G930100 |
| Latex free | yes | Preparation instructions | MPK02 |
|  | | | |
| Medicon media | | | |
| Complete picture | | Detail images | |
| | | | |
| size | Ø 1.45x18mm 27.06mm | Coating | |
| Storage conditions | | | |
| No special requirements | | | |
| Materials | | | |
| | Material number (DIN designation) | | |
| Main part | 1.4112 (X90CrMoV18) | | |