




## Product data sheet REF 68.79.05

11.12.2025

General data			
article number	68.79.05	Twist drill Ø1.45x6,5mm 155.6mm	
Unit	piece	Customs code (HS code)	90184910
GMDN	10520	GTIN	04046826540692
Risk category	2a (EU MDD)	UMDNS-Code	11-341
Classification	6 (EU MDD)	Base UDI-DI	
Net weight (kg)	0,000	Operating instructions	G930100
Latex free	yes	Preparation instructions	MPK02
			
Medicon media			
Complete picture		Detail images	
size	Ø 1.45x6,5mm 155,6mm	Coating	
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part	1.4112 (X90CrMoV18)		