



Product data sheet REF 68.80.30

11.12.2025

| General data | | | |
|-----------------|-------------|-----------------------------------|----------------|
| article number | 68.80.30 | titanium plate for eye prosthesis | |
| Unit | piece | Customs code (HS code) | 90213990 |
| GMDN | 46642 | GTIN | 04046826123604 |
| Risk category | 2b (EU MDD) | UMDNS-Code | 13-049 |
| Classification | 8 (EU MDD) | Base UDI-DI | |
| Net weight (kg) | 0,035 | Operating instructions | G688030 |
| Latex free | yes | Preparation instructions | MPK06 |



| | |
|---------------|-------------------------|
| Medicon media | 453.12.00 (Epiplating) |
|---------------|-------------------------|

Complete picture

Detail images



Storage conditions

No special requirements

Materials

| | Material number (DIN designation) |
|-----------|-----------------------------------|
| Main part | 3.7035 (ASTM_Grade_2) |