



Product data sheet REF 68.80.65

11.12.2025

General data			
article number	68.80.65	fastening screws p/2 f 68.80.64	
Unit	box/2	Customs code (HS code)	90213990
GMDN	47748	GTIN	04046826123826
Risk category	2b (EU MDD)	UMDNS-Code	16-055
Classification	8 (EU MDD)	Base UDI-DI	
Net weight (kg)	0,001	Operating instructions	G688030
Latex free	yes	Preparation instructions	MPK06
  			
Medicon media	453.12.00 (Epiplating)		
Complete picture		Detail images	
			
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part	3.7165 (ASTM_Grade_5_TiAl6V4_(ELI))		