











## Product data sheet REF 68.80.84

11.12.2025

| General data  |                                   |   |                |
|---|-----------------------------------|---|----------------|
| article number  | 68.80.84                          | magnetic insert Maxi Ø 5.1x1.6mm f prosthesis plate   |                |
| Unit  | piece                             | Customs code (HS code)  | 90213990       |
| GMDN  | 46642                             | GTIN  | 04046826387402 |
| Risk category   | N.A. (EU MDD)                     | UMDNS-Code  |                |
| Classification  | NA (EU MDD)                       | Base UDI-DI   |                |
| Net weight (kg)   | 0,000                             | Operating instructions  |                |
| Latex free  |                                   | Preparation instructions  | MPK06          |
|    |                                   |   |                |
| Medicon media   |                                   |   |                |
| Complete picture  |                                   | Detail images   |                |
|    |                                   |   |                |
| size  | Ø 5.1x1.6mm                       | Coating   |                |
| Storage conditions  |                                   |   |                |
|  15 °C - 30 °C   |                                   |   |                |
| Materials   |                                   |   |                |
|   | Material number (DIN designation) |   |                |
| Main part   | 3.7065 (ASTM_Grade_4)             |   |                |