



Product data sheet REF 68.82.57

11.12.2025


General data			
article number	68.82.57	Stiller elevator 3mm retrograde d/e 22cm	
Unit	piece	Customs code (HS code)	90189084
GMDN	35048	GTIN	04046826124274
Risk category	1r (EU MDR)	UMDNS-Code	15-218
Classification	6 (EU MDR)	Base UDI-DI	4046826;B;R004;05;aj;6A
Net weight (kg)	0,039	Operating instructions	GATDR00X
Latex free	yes	Preparation instructions	MPK03-1



Medicon media	450.00.03 (general catalogue the art of surgery) 450.00.62 (catalogue Medicon instruments Oral and Maxillofacial Surgery)
---------------	--

Complete picture

Detail images



size	3mm 22cm	Coating	
------	----------	---------	--

Storage conditions

No special requirements

Materials

	Material number (DIN designation)
Main part	1.4104 (X14CrMoS17)