








## Product data sheet REF 68.82.58

11.12.2025

General data			
article number	68.82.58	Stiller compressor instr Ø 4mm d/e 18cm	
Unit	piece	Customs code (HS code)	90189084
GMDN	35048	GTIN	04046826124281
Risk category	1r (EU MDR)	UMDNS-Code	15-218
Classification	6 (EU MDR)	Base UDI-DI	4046826;B;R005;04;ak;6M
Net weight (kg)	0,045	Operating instructions	GATDR00X
Latex free	yes	Preparation instructions	
 			
Medicon media	450.00.03 (general catalogue the art of surgery) 450.00.62 (catalogue Medicon instruments Oral and Maxillofacial Surgery)		
Complete picture		Detail images	
  			
size	Ø 4mm 18cm	Coating	
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part	1.4104 (X14CrMoS17)		