







## Product data sheet REF 68.83.03A

11.12.2025

General data			
article number	68.83.03A	Ti Micro slot screw 1.5x3.5mm grey self-tapp p/10	
Unit	box/10	Customs code (HS code)	90211090
GMDN	61837	GTIN	04046826344436
Risk category	2b (EU MDD)	UMDNS-Code	16-101
Classification	8 (EU MDD)	Base UDI-DI	
Net weight (kg)	0,010	Operating instructions	G687004
Latex free	yes	Preparation instructions	MPK06
  			
Medicon media	451.08.13 (brochure intraoral distraction devices f upper a lower jaws)		
Complete picture		Detail images	
			
size	1.5x3.5mm	Coating	
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part	3.7165 (ASTM_Grade_5_TiAl6V4_(ELI))		