







## Product data sheet REF 68.83.13A

11.12.2025

| General data  |                                     |  |                |
|---|-------------------------------------|--|----------------|
| article number  | 68.83.13A                           | Ti Micro slot screw 1.5x13mm grey self-tapp p/10 |                |
| Unit  | box/10                              | Customs code (HS code)                           | 90211090       |
| GMDN  | 61837                               | GTIN   | 04046826344634 |
| Risk category   | 2b (EU MDD)                         | UMDNS-Code                                       | 16-101         |
| Classification  | 8 (EU MDD)                          | Base UDI-DI                                      |                |
| Net weight (kg)   | 0,009                               | Operating instructions                           | G687004        |
| Latex free  | yes                                 | Preparation instructions                         | MPK06          |
|    |                                     |  |                |
| Medicon media   |                                     |  |                |
| Complete picture  |                                     | Detail images                                    |                |
|   |                                     |  |                |
| size  | 1.5x13mm                            | Coating  |                |
| Storage conditions  |                                     |  |                |
| No special requirements   |                                     |  |                |
| Materials   |                                     |  |                |
|   | Material number (DIN designation)   |  |                |
| Main part   | 3.7165 (ASTM_Grade_5_TiAl6V4_(ELI)) |  |                |