






## Product data sheet REF 68.83.20A

11.12.2025

General data			
article number	68.83.20A	Ti Micro slot screw 1.5x20mm grey self-tapp p/10	
Unit	box/10	Customs code (HS code)	90211090
GMDN	61837	GTIN	04046826344757
Risk category	2b (EU MDD)	UMDNS-Code	16-101
Classification	8 (EU MDD)	Base UDI-DI	
Net weight (kg)	0,009	Operating instructions	G687004
Latex free	yes	Preparation instructions	MPK06
  			
Medicon media			
Complete picture		Detail images	
size	1.5x20mm	Coating	
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part	3.7165 (ASTM_Grade_5_TiAl6V4_(ELI))		