







Product data sheet REF 68.83.23A

11.12.2025

General data			
article number	68.83.23A	Ti Micro slot screw 1.5x3mm self-tapp Dr Kluba P/5	
Unit	box/5	Customs code (HS code)	90211090
GMDN	61837	GTIN	04046826439514
Risk category	2b (EU MDD)	UMDNS-Code	16-101
Classification	8 (EU MDD)	Base UDI-DI	
Net weight (kg)	0,000	Operating instructions	G687004
Latex free	yes	Preparation instructions	MPK06
  			
Medicon media			
Complete picture		Detail images	
			
size	1.5x3mm	Coating	
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part	3.7165 (ASTM_Grade_5_TiAl6V4_(ELI))		